




**Cell repair kit for standard flow cell, for Agilent diode array and multiple wavelength detectors**



Cell repair kit for standard flow cell, for Agilent diode array and multiple wavelength detectors.

Miniatura	Nr-art.	Nazwa
	G1315-68712	Cell repair kit for standard flow cell, for Agilent diode array and multiple wavelength detectors

**Opis**

Cell repair kit for standard flow cell, for Agilent diode array and multiple wavelength detectors.

**Produkt dostępny do wyczerpania zapasów.**